

**AED Pro® Accessories
and Training Products**

ZOLL®



Electrodes



CPR-D-padz®

The revolutionary one-piece adult CPR-D-padz, which contains a pair of electrodes, ensures fast and accurate placement. It senses and reports the depth and rate of the chest compressions to the ZOLL AED Pro®. The CPR-D-padz electrodes have a 5-year shelf life—the longest of any disposable electrode.

Part #8900-0800



Pedi-padz® II

ZOLL's unique intelligent pediatric capability is built into the AED Pro. Pedi-padz II electrodes are designed for patients under the age of 8 or under 55 pounds. The AED can sense when pediatric electrodes are attached and will switch from an adult to a pediatric ECG analysis algorithm. The Pedi-padz II electrodes have a two-year shelf life.

Part #8900-0810



CPR Stat-padz®

Two-piece electrodes with Real CPR Help® capability. 24-month shelf life.

Part #8900-0400 (8 pack)

Part #8900-0402 (single pack)



ECG Electrodes

200-pouch case. (3 ECG electrodes/pouch)

Part #8900-0003



AED Pro ECG Cable AAMI

Lead II patient monitoring with 3-lead cable.

Part #8000-0838

Accessories



SurePower® Battery Pack

The SurePower Battery Pack is a rechargeable lithium-ion power source that is reliable for both infrequent users and full-time professional rescuers. With onboard self-testing, usage and charging history, and an actual "Runtime" charge indicator ("Runtime" — 30 minutes remaining, versus "Capacity" — 25% remaining), you can be completely sure that your AED Pro is adequately powered. Also compatible with the ZOLL E Series® Monitor/Defibrillator. Charger not included.

Part #8019-0535



SurePower Charging Station

Test, recondition, and charge four SurePower batteries at once with the SurePower Charging Station. This system maximizes the life-span and reliability of ZOLL rechargeable batteries. Compatible with the SurePower Battery Pack, and all previous ZOLL rechargeable batteries. Also compatible with the SurePower II Battery Pack (requires bay adapter, sold separately). Batteries not included.

Part #8050-0030

SurePower Single-Bay Battery Charger

Test, recondition, and charge a single ZOLL SurePower or SurePower II Battery Pack. Not compatible with NiCad or sealed lead-acid batteries. Batteries not included.

Part #8200-000100





**AED Pro Non-rechargeable
Lithium-ion Pack**

5-year maximum usable life.
Cannot be recharged.

Part #8000-0860



AED Pro Soft Case

Part #8000-0810



AED Pro Molded Vinyl Case

Provides spare battery compartment.

Part #8000-0832



Hard-shell Carry Case

Hard-shell case with foam cutouts.

Part #8000-0875

Training



ZOLL Travel Trainer

This streamlined portable trainer includes a flattened manikin silhouette, which features a built-in heart rhythm simulator and CPR-D-Demo pad for demonstrating the AED Pro. Easily folds into the small duffel bag provided for convenience when carrying.

Part #8008-0006



AED Plus Demo Kit

Equally ideal for demonstrating Real CPR Help® with the AED Pro, the AED Plus® Demo Kit comes with a CPR-D-padz Demo pad and a manikin in a large carry bag. (Clinical AED Pro and simulator not included)

Part #8000-0834

AED Pro Simulator

Simulates heart rhythm and 3-lead ECG for AED Pro demonstration using clinical unit.



Part #8000-0829



See-Thru CPR® Simulator

This simulator provides common arrhythmias and lead II ECG for training and demonstrations, with the added benefit of simulating See-Thru CPR®—the AED Pro's cutting-edge technology for reducing the duration of pauses during chest compressions. Also compatible with ZOLL's X Series® and R Series® defibrillators.

Part #8009-0751



CPR-D-Demo padz

Demonstration electrodes for use with clinical AED Pro or AED Plus (requires AED and simulator).

Part #8900-5007



Demo CPR Stat-padz

Demonstration electrodes for use with clinical AED Pro or AED Plus (requires AED and simulator).

Part #8900-0190



Demo Manikin

Head adjustable for airway management.

Part #8000-0835



ZOLL Medical Corporation, an Asahi Kasei Group company, develops and markets medical devices and software solutions that help advance emergency care and save lives, while increasing clinical and operational efficiencies. With products for defibrillation and monitoring, circulation and CPR feedback, data management, fluid resuscitation, therapeutic temperature management, and ventilation, ZOLL provides a comprehensive set of technologies that help clinicians, EMS and fire professionals, and lay rescuers treat victims needing resuscitation and acute critical care. For more information, visit www.zoll.com.

The Asahi Kasei Group is a diversified group of companies led by holding company Asahi Kasei Corp., with operations in the chemicals and fibers, homes and construction materials, electronics, and health care business sectors. Its health care operations include devices and systems for critical care, dialysis, therapeutic apheresis, transfusion, and manufacture of biotherapeutics, as well as pharmaceuticals, diagnostic reagents, and nutritional products. With more than 30,000 employees around the world, the Asahi Kasei Group serves customers in more than 100 countries. For more information, visit www.asahi-kasei.co.jp/asahi/en/.

**ADVANCING
RESUSCITATION.
TODAY.®**

**ZOLL Medical Corporation
Worldwide Headquarters**
269 Mill Road
Chelmsford, MA 01824
978-421-9655
800-804-4356
www.zoll.com

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.

©2017 ZOLL Medical Corporation. All rights reserved. Advancing Resuscitation. Today., AED Plus, AED Pro, CPR-D-padz, E Series, Pedi-padz, Real CPR Help, R Series, See-Thru CPR, SurePower, X Series and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries.

All other trademarks are the property of their respective owners.

MCN PP 1503 0124-05

ZOLL®