

# ELECTRODE PRODUCT GUIDE








Product Name	Differentiation Features	Location Used	Defibrillator
<b>CPR-D-padz</b>	 <ul style="list-style-type: none"> <li>Enables Real CPR Help® and See-Thru CPR®</li> <li>One-piece design simplifies proper positioning on the chest</li> <li>One size fits all, with ability to expand for larger chest</li> </ul>	<ul style="list-style-type: none"> <li>EMS &amp; Fire</li> <li>Public Access</li> </ul>	<ul style="list-style-type: none"> <li>M Series®</li> <li>E Series®</li> <li>R Series®</li> </ul> <ul style="list-style-type: none"> <li>AED Plus®</li> <li>AED Pro®</li> <li>X Series®</li> </ul>
<b>CPR Stat-padz</b>	 <ul style="list-style-type: none"> <li>Enables Real CPR Help and See-Thru CPR</li> <li>Multi-function, two-piece design ensures all therapies can be delivered</li> <li>Accommodates AP and AA pad placement</li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> <li>Med/Surg floors</li> <li>ER &amp; ICU</li> <li>EMS &amp; Fire</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul> <ul style="list-style-type: none"> <li>AED Plus</li> <li>AED Pro</li> <li>X Series</li> </ul>
<b>OneStep Basic</b>	 <ul style="list-style-type: none"> <li>Multi-function electrode</li> <li>Automated self-test with R Series</li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>OneStep Pacing</b>	 <ul style="list-style-type: none"> <li>OneStep Basic features plus:                             <ul style="list-style-type: none"> <li>OneStep Pacing</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>R Series</li> </ul>
<b>OneStep CPR (AP)</b>	 <ul style="list-style-type: none"> <li>OneStep Basic features plus:                             <ul style="list-style-type: none"> <li>Real CPR Help</li> <li>See-Thru CPR</li> <li>Anterior-Posterior placement</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>OneStep CPR (AA)</b>	 <ul style="list-style-type: none"> <li>OneStep Basic features plus:                             <ul style="list-style-type: none"> <li>Real CPR Help</li> <li>See-Thru CPR</li> <li>Anterior-Anterior placement</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>OneStep Complete</b>	 <ul style="list-style-type: none"> <li>All OneStep Basic, Pacing, and CPR capabilities</li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>R Series</li> </ul>
<b>OneStep Pediatric</b>	 <ul style="list-style-type: none"> <li>Designed for patients 0–8 years of age and less than 55 pounds (25 kg)</li> <li>Automatically adjusts joules setting to pediatric levels</li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>OneStep Pediatric CPR</b>	 <ul style="list-style-type: none"> <li>Designed for patients 0–8 years of age and less than 55 pounds (25 kg)</li> <li>Automatically adjusts joules setting to pediatric levels</li> <li>Real CPR Help</li> <li>See-Thru CPR</li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>Pedi-padz® II</b>	 <ul style="list-style-type: none"> <li>Connects directly to AEDs &amp; ALS devices</li> <li>Pediatric intelligence function with ZOLL AEDs</li> </ul>	<ul style="list-style-type: none"> <li>EMS &amp; Fire</li> <li>Public Access</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul> <ul style="list-style-type: none"> <li>AED Plus</li> <li>AED Pro</li> <li>X Series</li> </ul>
<b>Pedi-padz Multi-function Liquid or Solid Gel</b>	 <ul style="list-style-type: none"> <li>Multi-functional electrode for pediatric patients</li> <li>Available in liquid or solid gel</li> </ul>	<ul style="list-style-type: none"> <li>Hospital</li> <li>EMS &amp; Fire</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>Pedi-padz Radiolucent Solid Gel</b>	 <ul style="list-style-type: none"> <li>For use in cath and RP labs for procedures</li> <li>Virtually invisible to x-ray</li> </ul>	<ul style="list-style-type: none"> <li>Hospital                             <ul style="list-style-type: none"> <li>Cardiac cath lab</li> <li>Electrophysiology labs</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>Pro-padz® Multi-function Liquid or Solid Gel</b>	 <ul style="list-style-type: none"> <li>Available in liquid or solid gel</li> <li>Available with or without tray</li> <li>Used in cath and EP labs for procedures</li> </ul>	<ul style="list-style-type: none"> <li>Hospital                             <ul style="list-style-type: none"> <li>Cardiac catheterization lab</li> <li>Electrophysiology labs</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>Pro-padz Radiolucent Liquid or Solid Gel</b>	 <ul style="list-style-type: none"> <li>For use in cath and EP labs for procedures - Virtually invisible to x-ray</li> <li>Pre-Connect (solid only)</li> <li>Speed pack</li> </ul>	<ul style="list-style-type: none"> <li>Hospital                             <ul style="list-style-type: none"> <li>Cardiac catheterization lab</li> <li>Electrophysiology labs</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>Pro-padz Sterile</b>	 <ul style="list-style-type: none"> <li>For use in operating rooms</li> <li>Helps maintain sterile barrier</li> <li>Packaged in two pouches</li> <li>Available in 54" or 10-foot lead wires</li> </ul>	<ul style="list-style-type: none"> <li>Hospital                             <ul style="list-style-type: none"> <li>Operating rooms</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>Stat-padz®</b>	 <ul style="list-style-type: none"> <li>Multi-function</li> <li>Pre-Connect</li> <li>Speed pack</li> <li>Different pad shapes and colors</li> </ul>	<ul style="list-style-type: none"> <li>EMS &amp; Fire</li> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> <li>X Series</li> </ul>
<b>Stat-padz II</b>	 <ul style="list-style-type: none"> <li>Connects directly to AEDs and ALS devices</li> </ul>	<ul style="list-style-type: none"> <li>EMS &amp; Fire</li> <li>Hospital</li> </ul>	<ul style="list-style-type: none"> <li>M Series</li> <li>E Series</li> <li>R Series</li> </ul> <ul style="list-style-type: none"> <li>AED Plus</li> <li>AED Pro</li> <li>X Series</li> </ul>

See reverse side for proper cable/connector requirements.

# ELECTRODE GUIDE



	CPR-D-padz®	CPR Stat-padz®	OneStep™ Complete	OneStep™ CPR AP	OneStep™ CPR AA	OneStep™ Pediatric CPR
<b>AED Plus®</b>	 <ul style="list-style-type: none"> <li>Connects directly for Real CPR Help® and defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for Real CPR Help and defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>
<b>AED Pro®</b>	 <ul style="list-style-type: none"> <li>Connects directly for Real CPR Help, See-Thru CPR®, and defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for Real CPR Help, See-Thru CPR, and defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>
<b>X Series®</b>	 <ul style="list-style-type: none"> <li>Requires CPR Connector (PN# 8000-0370) for CPR feedback</li> <li>Connects directly for defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Requires CPR Connector (PN# 8000-0370) for CPR feedback</li> <li>Connects directly for defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for all capabilities</li> <li>No OneStep Pacing Capabilities</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for all capabilities</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for all capabilities</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for all capabilities</li> </ul>
<b>E Series®</b>	 <ul style="list-style-type: none"> <li>Requires CPR Connector (PN# 8000-0370) for CPR Feedback</li> <li>Connects directly for defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Requires CPR Connector (PN# 8000-0370) for CPR feedback</li> <li>Connects directly for defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>	<ul style="list-style-type: none"> <li>Not compatible</li> </ul>
<b>R Series®</b>	 <ul style="list-style-type: none"> <li>Requires R Series® CPR Connector (PN# 8009-0020) for CPR Feedback</li> <li>Connects directly for defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Requires R Series CPR Connector (PN# 8009-0020) for CPR Feedback</li> <li>Connects directly for defibrillation</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for all capabilities</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for all capabilities</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for all capabilities</li> </ul>	<ul style="list-style-type: none"> <li>Connects directly for all capabilities</li> </ul>
<b>Shelf Life</b>	5 yrs	2 yrs	2 yrs	2 yrs	2 yrs	2 yrs
<b>Quantity and Part Numbers</b>	Individual: 8900-0800-01	Individual: 8900-0402 8 per case: 8900-0400	Individual: 8900-0224-01 8 per case: 8900-0214-01	Individual: 8900-0223-01 8 per case: 8900-0213-01	Individual: 8900-0225-01 8 per case: 8900-0217-01	Individual: 8900-000219-01 8 per case: 8900-000220-01

(AP: anterior-posterior; AA: anterior-anterior)  
See reverse side for individual product details.